

In vitro Diagnosis of Drug Allergy: Current Status and Perspectives

Thomas A. Fleisher, M.D., FAACAI, FACAAI
National Institutes of Health
Bethesda, MD, USA

Dr. Fleisher has no conflicts of interest related to this presentation

Drugs as Immunogens

- Biologics: foreign macromolecules (*e.g.* antibodies, recombinant proteins) act directly as immunogen
- Drugs (non-biologics)
 - Hapten – drug (*e.g.* β -lactam antibiotics, quinidine) combines with a host macromolecule
 - Pro-hapten – processed drug (*e.g.* sulfonamides, phenytoin) combines with a host macromolecule
- Drugs can act directly to stimulate an immune receptor (pharmacologic interaction with immune receptors = p-i concept)

Use of *in vitro* Testing for Drug Allergy

- Testing in the setting of an immediate drug reaction
- Testing in the setting of a delayed drug reaction
- Testing on the horizon

Immediate Reaction to Drug

- Gell and Coombs type 1 reaction that occurs rapidly upon exposure to a specific drug
- Standard approach to evaluate is immediate skin testing (penicillin major and minor determinants are validated, other drugs ?)
- In vitro methods of evaluation include:
 - Tryptase to establish mast cell degranulation
 - Allergen (drug) specific IgE testing
 - Basophil activation test (BAT)

Tryptase Testing

- Mature tryptase reflects mast cell degranulation and is elevated in a systemic allergic reaction
- Current laboratory test most widely available measure total tryptase (not mature tryptase)
 - Released within 30-60 minutes following activation and half life is ~2 hours allows longer “testing window”
 - Levels above normal range (vary among labs: 10-11.4 ng/mL) are consistent with anaphylaxis (or increased mast cell numbers) but the sensitivity is not high
 - More sensitive test for anaphylaxis: mature tryptase level or a total tryptase rise over baseline of > 2 ng/mL


Allergen Specific IgE Testing

- *In vitro* “equivalent” of immediate skin testing
- Does not subject patient to risk and does not have a potential of inducing sensitization
- Limited range of drugs available impacts utility: β -lactams (penicilloyl G & V, ampicilloyl, amoxocilloyl), ACTH, cefator, ceftriazone, chlorhexidene, ethylene oxide, gelatin, insulin, neuromuscular blocking agents, tetanus toxoid)
- Tests generally have high specificity with lower sensitivity - **negative test does not rule out allergy**

Basophil Activation Test

- Test evaluates basophils present in either whole blood or separated mononuclear cells
- Validated for aeroallergens, hymenoptera venoms, foods, latex, some drugs (generally based on a generated drug-protein complex)
- Commercial assay (not FDA approved in USA): uses expression of CCR3 to identify basophils and expression of CD63 to identify activation after incubating cells with drug complex
- “Enhanced assay” adds a third marker, CD203c

Basophil Activation Test



Positive control: 52.5% CD63+, SI - $5501/386 = 14.2$

Positive drug BAT: 20.6% CD63+; SI - $1893/386 = 4.9$

Basophil Activation Test

- Advantages
 - Does not subject patient to any risks
 - Functional test that resembles the *in vivo* pathway
 - Relatively good sensitivity with high specificity
 - Positive BAT depends on type of allergen
 - Aeroallergens/foods $\geq 15\%$ CD63+ basophils
 - Venoms $\geq 10\%$ CD63+ basophils
 - Drugs (β -lactams, analgesics) $\geq 5\%$ CD63+ basophils
- Disadvantages
 - Must have viable, non-activated cells (24 hr “window”)
 - More limited availability since it requires a flow cytometer and generation of drug-protein (hapten-carrier) complex
 - **Negative test does not rule out drug allergy**

BAT in Radiocontrast Media Reactions

- Evaluation of 26 patients with history of immediate radiocontrast media (RCM) reactions: BAT using five different RCM products (tested months later)
- BAT results: 15/26 patients had a positive BAT
 - 1:100 RCM: patients = 13.1% CD63+/SI=8.1 (p=0.01)
controls = 2.7% CD63+/SI=1.5
 - 1:10 RCM: patients = 19.2% CD63+/SI=9.0 (p=0.001)
controls = 3.7% CD63+/SI=2.3
- Receiver Operator Curve (ROC) area under the curve was 0.79 = test with moderate accuracy

Delayed Immunologic Reaction to Drugs

- Most commonly linked to cellular response (Gell and Coombs Type IV reaction involving T cells)
- These reactions have been subdivided into
 - Type IVa: mediated by Th1 response
 - Type IVb: mediated by Th2 response
 - Type IVc: mediated by cytotoxic cell response
 - Type IVd: mediated by neutrophilic inflammation
- Additional data now suggests that some reactions involve conventional TcR activation (*e.g.* where there is an HLA link) and others involve direct drug-immune receptor interaction (p-i concept)


Focus of *in vitro* Testing

- Confirm that the clinical findings are the result of an immunologic response (rather than a pharmacologic or idiosyncratic response)
- Identify the causative drug in settings where multiple drugs have been administered
- Current testing methods
 - Lymphocyte transformation test (LTT)
 - CD69 upregulation flow cytometry test
 - Cytokine production
 - Evaluation of cytotoxicity (or its products)

Lymphocyte Transformation Test (LTT)

I- Activation *in vitro*

Varied concentrations of pure drug, incubate at 37°C with 5% CO₂



II- Quantify Response

Harvest cells and count radioactivity,
results: cpm or stimulation index (SI =
drug stimulated cpm/unstimulated cpm)

Lymphocyte Transformation Test (LTT)

- Must use controls to establish lack of drug induced toxicity and to rule out non-specific activation
- Must have viable cells and requires sterile tissue culture
- LTT has been successfully applied to drug associated:
 - Maculopapular exanthem
 - Pustular exanthem
 - Stevens Johnson syndrome (SJS)/toxic epidermal necrolysis (TEN)
 - Drug rash with eosinophilia and systemic symptoms (DRESS)
- Positive LTT has generally been defined as a stimulation index ($SI = \text{cpm with drug} / \text{cpm with medium}$) > 2
- Sensitivity is 60-70% under optimal conditions with a higher specificity
- **Negative test does not rule out T cell mediated drug response**

Evaluation of LTT in Different Types of Delayed Hypersensitivity Drug Reactions

- 27 patients in three groups: 8 maculopapular eruptions (MP), 6 SJS + 2 TEN, 11 DRESS
- Evaluated by LTT at 1 week, 2-4 weeks, 5-8 weeks, 1 year and > 1 year following onset
- Patients with MP and SJS/TEN had positive LTT at 1 week post-onset, response declined over time
- Patients with DRESS were negative at 1 week and were positive at 5-8 weeks

LTT Used to Identify the Drug that Induced DRESS

- Two patients receiving multiple drugs including anticonvulsants and antibiotics associated with the development of DRESS
- Evaluation by LTT utilized all drugs that had been given, each at 7 concentrations (1-200 $\mu\text{g/ml}$)
- Studied 3 months after the clinical presentation
- Causative drug was identified as ceftriaxone in one pt and piperacillin-tazobactam in the other pt
- LTT assay proved valuable in defining the drug associated with DRESS (avoid in the future)

LTT Summary

- LTT appears to be a suitable complement to other testing in delayed drug reactions
- Time line of positivity may differ between the different types of delayed drug reactions
- Positive test helps identify the offending drug but a **negative test does not rule out drug related hypersensitivity**
- The test remains a research tool, it is not standardized and it requires tissue culture with results available after six or more days


Alternatives to LTT (^3H Thymidine)

- Evaluation of upregulation of a T cell activation antigen in response to *in vitro* drug exposure
 - CD69 up-regulation, an early product of T cell activation, measured by flow cytometry at 48 hrs
- Ex vivo cytokine production
 - Cytokine secretion into the supernatant following mononuclear cell culture with drug (*e.g.* γ -IFN)
 - Elispot assay measures individual T cell production of a cytokine following *in vitro* drug stimulation

T cell CD69 Upregulation


I- Activation *in vitro*

Varied concentrations of pure drug, incubate at 37°C with 5% CO₂



Evaluate T cells by
flow cytometry

II- Quantify Response




CD69 upregulation expressed as
percent CD69 positive T cells

CD69 Upregulation in Response to Drug

Evaluation of a
phenytoin-allergic patient
following 48 hrs of
stimulation
medium - negative control
Tetanus toxoid - positive control

Phenytoin - positive test

Unrelated drug clonazepam -
negative test




Summary of LTT Alternatives

- CD69 upregulation appears to perform similar to LTT with the advantage of being a 48 hour assay and not requiring radionuclides
- Cytokine production assays correspond to LTT but the actual cytokine produced does not appear to correlate well with the clinical phenotype (*i.e.* IFN- γ is typically produced with all types of delayed drug reactions)

Immunopathogenesis of SJS/TEN

Bullous skin processes (SJS/TEN) associated with drugs appear to be linked to cytotoxic T cell activity



Soluble Fas ligand (sFasL) and granulysin have been found in the serum of patients with SJS/TEN

“Real Time” Test to Diagnose SJS/TEN

- The serum level of granulysin is ~100X greater than sFasL in SJS/TEN making it an attractive target
- An immunochromagraphic test for serum granulysin (≥ 10 ng/mL) predicted SJS/TEN 2-4 days prior to mucocutaneous reuptions
- This assay could prove useful in predicting when a drug reaction will lead to SJS/TEN

In the Future

- Multiplex cytokine evaluation following in vitro culture (e.g. IFN- γ , IL-2, IL-4, IL-5, IL-8, IL-13, IL-17, etc) may reveal specifics about the type of immune response
- Nature of drug derived epitopes inducing an immune reaction often are not well understood
 - Mass spectrometry (MS) has evolved as a powerful tool to evaluate proteomics and metabolomics
 - MS used to characterize the functional antigens derived from piperacillin (in CF patient serum) with the identification of multiple drug derived haptenic structures bound to albumin
(Whitaker P, et al. J Immunol 2011, 187:200)

Summary *in vitro* Testing in Drug Allergy

- Immediate drug reactions
 - Specific IgE testing: safe test but there are limited numbers of suitable drug conjugates available for testing
 - BAT: promising functional test that requires viable cells and a drug conjugate preparation for activation
- Delayed drug reactions
 - Lymphocyte transformation test (LTT)
 - Most common research method to determine responsible drug
 - Issues remaining include: standardization, requirement for viable cells, six day sterile tissue culture period and use of radionuclides
 - CD69 upregulation may be equivalent to LTT – under study
 - *In vitro* cytokine production to drug – under study
 - Product of cytotoxic cells (granulysin) promising to help dx SJS/TEN prior to mucocutaneous symptoms (further study)

Conclusions

- The clinical story remains the most important starting point evaluating possible drug allergy
- *In vitro* testing can be complementary to *in vivo* testing and is evolving for the evaluation of both immediate and delayed drug allergy
- There is currently no single laboratory test that reliably establishes the drug responsible for an immunologically mediated drug reaction

References

- Fujita Y, et al. Rapid immunchromatographic test for serum granulysin is useful for the prediction of SJS/TEN. *J Am Acad Dermatol*. 2011, 65:65.
- Hausmann OV, et al. The basophil activation test in immediate-type drug allergy. *Immunol Allergy Clin North Am*. 2009, 29:555.
- Kano Y, et al. Utility of the LTT in the diagnosis of drug sensitivity: dependence on its timing and the type of drug eruption. *Allergy*. 2007, 62:1439.
- Lochmatter P, et al. In vitro tests in drug hypersensitivity diagnosis. *Immunol Allergy Clin North Am*. 2009, 29:537.
- Pichler WJ, et al. Immune pathogenesis of drug hypersensitivity reactions. *J Allergy Clin Immunol*. 2011, 127:S74.
- Romano A, et al. Diagnosis and management of drug hyper-sensitivity reactions. *J Allergy Clin Immunol*. 2011, 127:S67.